


# VITA MFT

Multi Functional Teeth

## Sales Information




VITA shade taking

VITA shade communication

VITA shade reproduction

VITA shade control

Date of issue 01.09.




**VITA**

Exclusively available in the  
VITA SYSTEM 3D-MASTER shades

## VITA MFT The perfect tooth for all circumstances

Multi Functional Teeth

### New standards with VITA MFT



Bar code and CE 0124 mark on the back

- Suitable for every user – simple, efficient and predictable set-up procedure
- Universal usage – one tooth for every indication
- Small shade/mould range for high economic efficiency
- VITA SYSTEM 3D-MASTER including “Best Match VITAPAN classical” for secure shade selection

### Competitive edge



- Attractive design for dentures staying true to the original
- Lifelike layering for a vivid play of colors
- Densely cross-linked and free of residual monomers which leads to outstanding biocompatibility
- High crack and abrasion resistance for long durability




### Functional principle

#### **Convex design of labial facettes for a systematic lip contour**

The slightly convex anterior teeth support and promote dentofacial harmony and life-like aesthetics.

#### **Slightly broadened tooth neck areas with anteriors and posteriors**

Slightly broadened cervical areas at anteriors and posteriors reduce the effort for a true-to-nature and age-appropriate design of the interdentium.



#### **Design of palatal surfaces for the normal phonetic behavior**

Through the nature-identical shape of the anteriors the patient is able to maintain his/her phonetic skills and the tongue contact he/she is used to.

#### **Reduced ridge lap for less grinding and more bonding surface**

Reduced ridge lap allows for a significant reduction of grinding work and enlarges the surface for safe bonding to the denture base material.



#### **easy-centric® for easy occlusion determination**









The perfect harmony of the tooth design ensures professional as well as economical set-up procedures for beginners and pros alike.

#### **Pre-abraded occlusal surfaces for minimal shear forces**

The specifically designed cusps support an ideal static alignment of the prosthesis. This is particularly important with implant-supported prostheses.



## Shades

VITA SYSTEM 3D-MASTER							
							
1M1	2L1.5	2M1	2M2	3L2.5	3M2	3R2.5	4M2
(~ B1)*	(~ B2)*	(~ A1)*	(~ A2)*	(~ B3)*	(~ A3)*	(~ A3.5)*	(~ A4)*

\* Best Match VITAPAN classical

## Selection

- 4 upper anterior sets – Triangular
- 3 upper anterior sets – Ovoid
- 3 upper anterior sets – Rectangular
- 2 upper anterior sets – Square
- 5 lower anterior sets (1 Long Neck)
- 3 upper posterior sets
- 3 lower posterior sets



## Indications

- Full dentures
- Implant prosthetics
- Telescope and cone constructions
- Attachment Dentures
- Model cast prostheses
- Interim solutions

# VITA MFT Convincing performance all along the line

Multi Functional Teeth



Wall/tabletop mould chart, Interactive mould guide, Packing, Product information, User manual, Mould chart, Ads

## Product information

- Essential information on VITA MFT summarized for the customer

## User Manual

- Interesting facts about prosthetics in general and about the usage of VITA MFT in detail: a meaningful guide for users and beginners

## Mould chart

- Overview of the VITA MFT shade and mould range including combination table and available forms

## Wall mould chart

- VITA MFT at a glance


## Interactive mould guide

- All VITA MFT sets at a glance – ideal for the tooth selection directly with the patient

## Ads

- Instruments for effective press releases regarding VITA MFT

With the unique VITA SYSTEM 3D-MASTER, all natural tooth shades can be systematically determined and perfectly reproduced.



Please note: Our products were developed specifically for dental use and must be processed according to the instructions for use. The information/illustrations contained herein do not constitute an assurance of properties and are not binding.

VITA Zahnfabrik is certified according to the Medical Device Directive; the following products bear the mark  0124 : **VITA MFT**

US 5498157 A · AU 659964 B2 · EP 0591958 B1

# VITA

VITA Zahnfabrik H. Rauter GmbH & Co.KG  
Postfach 1338 · D-79704 Bad Säckingen · Germany  
Tel. +49(0)7761/562-0 · Fax +49(0)7761/562-299  
Hotline: Tel. +49(0)7761/562-222 · Fax +49(0)7761/562-446  
[www.vita-zahnfabrik.com](http://www.vita-zahnfabrik.com) · [info@vita-zahnfabrik.com](mailto:info@vita-zahnfabrik.com)